GB Dublisil

Addition-vulcanising vinyl-polysiloxane in several final hardnesses for all dental duplications. Easy and precise to mix in a ratio of 1:1.

Product description:

Dublisil products are universal duplication materials for all kinds of duplications. Due to the practical way of mixing and the very thin-flowing initial consistency absolutely price duplication models can be produced. Dublisil forms are nearly shrinkless and have a high form consistency as well as a long shelf life due to the used platinum hardener system. Negative form made of Dublisil can easily and poured out without losing precision for several times. The surface of the duplication model is smooth and hard without after-treatment.

Preparatory measures:

Almost every material e. g. plaster models, model resins and metals can be duplicated with Dublisil.

Important: The master model must not be watered but should be dry and free of any insulating agents etc. For an economical consumption of the duplication material we recommend to make use of the Dublisil flask (item-no. 3805) with a variably adjustable interior as well the foam space maintainer (item-no. 38076) for the blocking out of the base parts which are unimportant for the duplication.

Dosage:

Dublisil is mixed in the ratio 1:1 of the components A (catalyst/white) and B (base/coloured).

Processing:

The two single components can manually be mixed in a mixing cup or automatically by means of a dosage and mixing unit. The most rational and clean method is the mixing in the Dosper M/4. At the touch of the button Dublisil is fully automatic and losslessly measured out, mixed and produced without air bubbles. In case that the material is manually mixed the fill in the material in a preferably thin stream into the flask.

When using the Dosper M/4 this stream is defined by the unit. The distance to the flask should be approx. 20 cm.

Further processing:

The best results regarding precision and homogeneity are achieved by the so-called pressure duplication technique. For this the just filled flask is put into a pressure pot (e. g. Pneuress T or Polymax). In case that the duplication form is produced under pressure it is important that also the model is produced under the same pressure. After the vulcanisation has been completed the master model can directly be deflasked. Before pouring out the duplication form is moistened with Spannex II. This leads to a lower surface tension and improves the flow properties of the model material. All dental plasters, model resins and embedding materials are suitable.

General aspects:

Dublisil is unproblematic in physiological and toxicological respects. A chemical and thermal disinfection is possible. Wear protective clothing! Silicone impurities on clothes cannot be cleaned - also not chemically. Allergic reaction to product components cannot be excluded in case of sensitive persons. The vulcanisation can be disturbed by latex gloves, hand cremes, cleaning agents etc. The product is exclusively designed for the a. m. field of application and is only to be used by qualified and instructed persons.

The general information of this product is based on detailed application practice and intensive research work. Due to further development we reserve the right of modification.

Technical data	Dublisil 15	Dublisil 20	Dublisil 30
Mixing time	45 sec. ± 15 sec.		
Processing time	6 ± 1 min.		
Deflasking time	25 ± 5 min.		
Strain in compression	19.2 %	14.3 %	11.0 %
Recovery from deformation	99.9 %	99.8 %	99.8 %
Linear dimensional change	<0.1 %		
Elongation	305 %	235 %	180 %
Colour	green	pink	lilac
Final hardness	15 shore A	20 shore A	30 shore A

The mixing and processing time refer to a room temperature of 23 °C (73.4 °F) and a relative air humidity of 50 %.

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