Uso previsto

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Intended use

Dental thermoforming blank for the production of bandage plates, drill guides and various splints.

Field of application

Biolit D is especially suitable for the production of therapy splint – added with resin – and for the production of drill guides in the implantation technique.

IMPORTANT: Due to the special chemical composition and the corresponding very inert temperature absorption **Biolit D** can be sterilised without problems. **Biolit D** thermoforming foils can be individualised excellently with all PMMA resins as well as when using a bonding with a light-curing resin.

Preparation

Biolit D on the based on polycarbonate has to be predried before use. The air humidity which is absorbed during the storage partly causes considerably blister forming without predrying. In order to receive results without any blisters **Biolit D** has to be predried for approx. 4 hours at a temperature of 80 °C (176 °F) per each mm layer strength. The permanent storage in the drying chamber is not recommendable as the materials react with a slight yellow discoloration.

Processing

Please take the heating times of Dreve thermoforming units from the following scale:

Biolit D	1.0 mm	2.0 mm
Drufomat	1:45 min	2:15 min
Drufosmart	1:45 min	2:15 min
Drufomat SQ / TE	1:45 min	2:15 min
Vacfomat U / P	1:45 min	2:15 min
Vacfomat 2000	1:35 min	2:05 min

Before thermoforming the plaster models should be watered and insulated. Through vacuum or effect of pressure the soft plasticised blank is formed. The cooling down phase under pressure / vacuum should correspond at least to the particular heating time. **Biolit D** can be reworked with carbide cutters or scissors. For polishing of the cut special flexible grinding disks (**UltraTrim**) are suitable.

Delivery forms

Ø 120 / 125 mm	clear transparent
□ 120 / 125 mm	clear transparent
Refractory deformation temperature:	150–160 °C (302–320 °F)
Pyrolysis:	> 230 °C
Density:	1.2 g/cm ³

REF D420013	1.0 mm	10 pcs.
REF D420014	2.0 mm	10 pcs.

This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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