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Dental use
Intended splinting material short-term

Extended purpose
Dental thermoforming blank for the production of moulded parts for the crown and bridge technique.

Indications for use
Drufolen H is especially suitable for the production of hard-flexible cast form parts for the C + B technique and for the template technique at C + B temporaries. Drufolen H does not bond with acrylate and burns without residues.

Contraindication
In case of allergies to any of the ingredients do not use Drufolen H. Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e. g. allergies) or irritations principally cannot be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

Target patient group
Patients who need dental care

Intended users
Trained staff in dental field

Technical data
Technical datasheet available on request.

Processing
Before the thermo-forming process the plaster models should be watered and insulated with an insulating agent on alginate basis (Isolat film). By means of vacuum or pressure the plastified plate is formed on the object. The cooling down phase under pressure / vacuum should at least correspond to the respective heating time. Drufolen H can be reworked by means of a scissors. For the smoothing of the incisions special flexible grinding wheels (UltraTrim) are suitable. The processing with a hand adapter is also possible. For the production of templates for the c & b temporaries the blank is thermo-formed on the situation model and cut

The product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self

in the areas of the of the ungrounded adjacent teeth which are located in the area which should be cared.

After the preparation the blank template can be repositioned in the patient’s mouth resp. on the plaster model. The exact fixation is effected above the remaining teeth. The hollow space which results in the area of the preparation is filled with a temporary c & b resin (e. g. Fixtemp) and put back into the mouth / onto the insulated model. Please refer to the instruction of use of the corresponding resin materials for any further information regarding the processing.

Desinfection
A disinfection with all commercial materials is possible.

Note
• While working out on the material, it is recommended to use an exhaust system and protective glass or goggles as well as a face mask

Serious incidents
All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/ or patient is settled.

Lot number / Expiry date
The lot number and the expiry date are indicated on each Drufolen W packaging. In case of claims, please always indicate the lot number of the product. Do not use the product after the expiry date.

Storage conditions
At room temperature, i. e. 18°C–28°C.

Disposal
Disposal of contents/container in accordance with local/ regional/national/international regulations.

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denté produit avant using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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