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Material description

Super-transparent, slightly flexible thermoforming blank made of polyethylene for the production of dental splints. Bonds with acrylate.

Field of application

Suitable for the production of stable bleaching and medication splints ending at the gingival margin and for situations with slight undercut areas.

Processing

BioBleach® hard will be delivered with an insulating foil. This extremely thin foil (0.07 mm) is thermoformed together with the **BioBleach® hard** on the model side and guarantees a decisively improved transparency of the splint. Thus an additional insulation of the model surface is no longer necessary.

If processed in a Vacfomat, the insulation foil must be removed before thermoforming.

The adjustment of the heating and cooling times can either be taken from the survey of thermoforming foils, or can be programmed directly into the unit with the help of the barcode on the **BioBleach® hard** packaging.

Before thermo-forming the blank a space maintainer as a reservoir for the bleaching material has to be applied on the labial areas of the teeth.

a) With the light curing block out material **Cenit L** (REF D3815) the areas of the plaster model, which are supposed to be bleached, are evenly covered and afterwards get cured.

b) The whole model is covered with the thermo-forming blank **Drufolen W** 0.8 mm (REF D423481N). Then the part of the model to be bleached (incl. a slight overlapping of the cutting edges) is cut out of the blank by means of a pair of scissors and put back onto the model again.

IMPORTANT: In order to avoid any contact between the mucosa and the bleaching material the space maintainer should end approx. 1 mm above the marginal gingiva.

Keep away from children! Only for dental use!

Storage conditions

At room temperature, i. e. 18–28 °C.

Durability: see date of expiry.

BioBleach® hard is not to be used when expiry date has passed.

Protect against humidity, always close foil bag when packs are in use.

Delivery form

Ø 120 mm, clear transparent (special measures can be delivered on request).

REF	D420015	BioBleach® hard
		0,5 mm, 25 pcs.

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This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible

use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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