

Phrozen Sonic CS+

Qualified by Dreve



For all, who simply want to print

Compact LCD desktop printer with 22 µm resolution and an innovative 385 nm LED light source, delivering high light intensity for fast build times. The high-quality full-metal housing and improved vat and platform handling make the Phrozen Sonic CS+ a user-friendly 3D printer for all dental applications.

Features

- + High light intensity for short building times
- + Print process monitored via webcam
- + Computer-guided full control over the printer
- + Compatible with the entire FotoDent® product range
- + Easy replacement of wear parts
- + Qualified by Dreve: Defined settings for optimal print results
- + Maintenance and repair at Dreve locations
- + Inclusive: Building platform, vats, USB-Stick, spatula, funnel, and Slicing Software Dreve Elements



Highly economic
Lowest price




Highly compatible
Optimally matched to the digital workflow



Super simple
User-friendly handling and comprehensive support from experts



Super fast
Individual components are printed in less than 2 hours



Type	LCD Printer
Wavelength	385 nm
Light intensity	5 mW/cm ²
Auflösung	22 µm
Print volume	Standard: 74 x 168 x 185 mm ³ Mini: 75 x 75 x 185 mm ³
Vat volume	Standard: 350 ml, Mini: 150 ml
Dimensions	550 x 380 x 360 mm ³ (H x W x D)
Layer thickness	0.01–0.30 mm